



DTi FabTaq HS

Instruction Manual for use



User Manual for DTi FabTaq HS

Catalogue number: DT0102.250

Note: Applicable to all pack sizes

Manufactured by
DSS Takara Bio India Pvt Ltd
A-5 Mohan Co-operative, Industrial Estate,
Mathura Road, New Delhi, Delhi 110044

1) Product Description:

DTI FabTaq HS is a hot start PCR enzyme derived from FabTaq that includes a neutralizing monoclonal antibody that recognizes Taq DNA polymerase. This antibody binds to Taq polymerase and prevents non-specific amplification due to mispriming and/or formation of primer dimers before thermal cycling. The antibody is denatured during the initial DNA-denaturation step, allowing this product to be used with standard PCR conditions.

2) Components:

Components	Cat# DT0102.250 (250 U)
DTI FabTaq	250 Units
10X PCR Buffer (Mg ²⁺ plus)	1 ml
dNTP Mixture	800 µl

3) Concentration: 5U/ul**4) Storage Buffer:**

Tris-HCl (pH8.0)	20 mM
KCl	100 mM
EDTA	0.1 mM
DTT	1 mM
Tween 20	0.5%
NP-40	0.5%
Glycerol	50%

5) Supplied 10X PCR Buffer (Mg²⁺ plus)

Composition (10X)

Tris-HCl (pH8.9)	100 mM
KCl	500 mM
MgCl ₂	15 mM

6) Supplied dNTP Mixture (2.5 mM each):

dNTP Mixture is ready for use in PCR without dilution.

Form : Dissolved in water (sodium salts), pH 7 - 9

Purity : \geq 98% for each dNTP.

7) Storage: -20°C**8) Unit definition:**

One unit is the amount of enzyme that will incorporate 10 nmol of dNTPs into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

9) Purity:

Nicking, endonuclease, and exonuclease activity were not detected after incubation of 0.6 µg of supercoiled pBR322 DNA, 0.6 µg of λDNA, or 0.6 µg of λ-Hin d III digest with 10 U of this enzyme for 1 hour at 74°C.

10) Applications:

- For DNA amplification by hot start PCR
- For DNA sequencing

11) PCR products:

As most PCR products amplified with DTI FabTaq HS have one A at the 3'-termini, the obtained PCR products can be directly cloned into a T-vector. It is also possible to clone the product in blunt-end vectors after blunting and phosphorylation of the ends.

12) Materials required but not provided

Reagents: PCR primers, sterile purified water, template

Equipment: Thermal cycler (DTI FabSpeed thermal cycler)

Consumables: PCR tubes, Micropipettes and tips

13) Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot.

14) Protocol

General reaction mixture for PCR (total 50 µl):

Reagent	Volume
DTI FabTaq HS (5 U/µl)	0.25 µl
10×PCR Buffer (Mg ²⁺ plus)	5 µl
dNTP Mixture (2.5 mM each)	4 µl
Template	< 500 ng
Primer 1	10 - 50 pmol (final conc. 0.2 - 1.0 µM)
Primer 2	10 - 50 pmol (final conc. 0.2 - 1.0 µM)
Sterile purified water	up to 50 µl
Total	50 µl

NOTE: Reaction mixtures can be set up at room temperature. Be sure to keep all reagents on ice.

15) PCR Conditions:

This enzyme can be used with standard PCR conditions, since the monoclonal antibody is denatured in the initial DNA-denaturation step. There is no need for an additional step to denature the anti-Taq antibody.


Example: Amplification of a 1 kb DNA fragment

98°C	10 sec	} 30 cycles
55°C*	30 sec	
72°C	1 min	

NOTE: Denaturation conditions vary depending on the thermal cycler and tubes used for PCR. Denaturation for 5 - 10 sec at 98°C, or 20 - 30 sec at 94°C is recommended.

16) Experimental sample:

a) Amplification data: Good performance of this product has been confirmed by PCR using λ DNA as template (amplified fragment size 8 kbp). The protocol used for the assay and the results are as follows:


**b) Hot start function data:**

Test target sample added to a PCR system containing FabTaq and FabTaq HS respectively. Following an incubation at RT for 30 mins, the products were amplified using a thermal cycler. Non-specific bands (280 bp) were observed in the PCR system along with specific target bands (838 bp) using FabTaq, while there were no non-specific bands observed in PCR system using FabTaq HS.

The protocol used for the assay and the results are as follows:

Incubate at 25C for 30 mins
↓
Place the tubes in thermal cycler and run the program

94C	30 sec	} 25Cycles
55C	30 sec	
72C	30 sec	
72C	5 min	



Lane M: 100 bp DNA marker
Lane 1-3: FabTaq PCR system
Lane 4-6: FabTaq HS PCR system

Visit <https://store.dsstakarabio.com/pages/dti-fabtaq-hs> for more detailed product information






For more information contact directly below;

Address: DSS Takara Bio India Pvt Ltd, A-5 Mohan Co-operative, Industrial Estate, Mathura Road,
New Delhi, Delhi 110044

Email: enquiries@dsstakarabio.com

Toll-Free number 1800-212-4922

Description of Symbol Used:

-  Catalogue number
-  Batch Code
-  Date of Manufacturing
-  Use-by-date
-  Contains sufficient for <n> tests