

DTI FabTaq & DTI FabTaq HS

DTI FabTaq

DTI FabTaq DNA Polymerase is a recombinant version Taq polymerase derived from the *Thermus aquaticus* YT-1 strain, and is suitable for routine PCR applications.

For individual reaction setup and optimization, use individual components of enzyme, dNTP mixture, and 10X reaction buffer (with or without Mg^{2+}). A convenient 2X PCR master mix of Takara Taq enzyme, buffer, and dNTP mixture is available as Premix Taq DNA Polymerase.

- Excellent for standard PCR amplifications.
- Low bacterial DNA contamination.
- Premix Taq enables easy PCR setup and minimizes pipetting steps




Made in India

Application

- Routine PCR
- DNA sequencing

94C 1 min
94C 30 sec
65C 10 min
4C ∞ } 30 cycles



Lane M: λ Hind III digest marker
Lane 1-3: Test in triplicates

Figure 1: Gel electrophoresis image of PCR

Amplification consistency of 8 kb DNA fragment
Consistent amplification of 8 kb fragment from λ DNA is observed using the DTI FabTaq. The

Want to know more?
Scan QR code



Why Choose DTI Brand ?

DTI is our new in-house brand, we are manufacturing products in India, where we can design products for you and deliver them to you at affordable prices.



ISO 9001 certification



Quick Turn around time



Affordable Price



Strict temperature
Control



Product customization

DTI FabTaq HS

An antibody-mediated hot-start version of DTI FabTaq Polymerase, which is a recombinant version of full-length Taq polymerase. It has the same characteristics and capabilities as the native Taq polymerase, and is suitable for a variety of standard PCR applications. The polymerase is supplied with separate tubes of buffer (Mg²⁺ plus) and dNTPs.

Application

- Hot-start PCR
- Multiplex PCR

Comparison with a standard Taq DNA polymerase


The performance of DTI FabTaq HS in a multiplex reaction was compared to its performance in single-target reactions, and to that of the standard DTI Fab Taq DNA Polymerase in a multiplex reaction (Figure 1). DTI FabTaq and DTI FabTaq HS were each used to amplify a human genomic DNA template with eight different primer pairs, each specific for a target ranging from 84 to 432 bp in size. In addition, Takara HS Taq was used to perform an individual amplification reaction with each of the eight primer pairs.

The sample amplified using DTI FabTaq HS (Lane 10) shows more efficient multiplex amplification of all the bands individually amplified in Lanes 1–8 than the sample amplified using DTI FabTaq (Lane 9), and it lacks the low molecular weight non-specific amplification band seen with DTI FabTaq (Lane 9). The intensities of the bands in Lane 10 are comparable to the intensities of the corresponding individual bands in Lanes 1–8, indicating that DTI FabTaq HS provides multiplex efficiencies comparable to those observed in individual amplification reactions.

Want to know more?
Scan QR code



Comparing the ability of DTI FabTaq HS to amplify human genomic DNA fragments to that of the standard Taq DNA Polymerase. PCR reactions were performed using human genomic DNA as a template and primer pairs for eight different targets. Lanes 1–8 contain individual reactions for each primer pair amplified using DTI FabTaq HS. Lanes 9 and 10 contain multiplex PCR reactions performed with all eight primer pairs in a single tube, amplified with either the standard Taq DNA polymerase (Takara Taq, Lane 9) or DTI FabTaq HS (Lane 10).



Product	Cat. No.	Pack size
DTI FabTaq	DT0101.250	250 Reactions
	DT0101.500	500 Reactions
	DT0101.1K	1000 Reactions
DTI FabTaq HS	DT0102.250	250 Reactions

DSS Takara Bio India Pvt. Ltd.

A-5, Mohan Cooperative Industrial Estate, Mathura Road, New Delhi 110044, INDIA

Toll Free: 1800-212-4922

E: enquiries@dsstakarabio.com W: www.dsstakarabio.com