




Get the clarity with colored protein ladders



DTI Low Range Colored Protein Ladder



DTI Broad Range Coloured Protein Ladder



Product	Cat. No.	Pack size
DTI Low Range Colored Protein Ladder	DT2501.100	100 lanes
DTI Broad Range Colored Protein Ladder	DT2502.100	100 lanes

DTI Protein Ladder are three-color protein standard as pre-stained proteins. Low Range OG 10 Colored Protein Ladder - molecular weights from 10 to 180 kDa. Broad Range OG 12 Colored Protein Ladder - molecular weights from 10 to 245 kDa. Proteins are covalently coupled with a blue chromophore when separated on SDS-PAGE (Tris-Glycine Buffer). The ladders are supplied in the gel loading buffer and is ready to use.



Why Choose DTI Brand ?

DTI is our new in-house brand, we are manufacturing products in India, where we can design products for you and deliver them to you at affordable prices.



ISO 9001 certification



Quick Turn around time



Affordable Price



Strict temperature control



Product customization

Want to know more?
Scan QR code



DSS Takara Bio India Pvt. Ltd.

A-5, Mohan Cooperative Industrial Estate, Mathura Road, New Delhi 110044, INDIA

Toll Free: 1800-212-4922

E: enquiries@dsstakarabio.com W: www.dsstakarabio.com